

New  
Innovation

# Drummond

## Portable Pipet-Aid® XP2

### The Evolution of Xcellence Continues

- New Ergonomic Design—The most comfortable pipettor you ever laid a hand on
- New Power Source for Uninterrupted Extended Operation—Can be charged while in use
- New Ultra Quiet Precision Pump—Great control for aspiration or dispensing

The new Drummond Pipet-Aid XP2 advances the industry standards for pipetting **control, comfort, and convenience**. Building on the years of experience of providing state-of-the-art tools that make your job in the laboratory safer, easier, and more convenient, the Pipet-Aid XP2 delivers more. The ergonomic design is more comfortable to handle. The new power source enables more convenient extended operation and the new more powerful pump provides great control and quieter operations. The unit is supplied complete with a power supply/charger, 4 extra filters, and a removable stand.

### DRUMMOND PORTABLE PIPET-AID® XP2

Cat. No.	Description
4-000-501	Pipet-Aid XP2, 110V
4-000-501-A	Pipet-Aid XP2, 240V, CE, Australia
4-000-501-I	Pipet-Aid XP2, 230V, CE, International
4-000-501-UK	Pipet-Aid XP2, 240V, CE, UK

### Replacement Parts

Cat. No.	Description
4-000-520	Nosepiece Cover—Front Cover—Filter Access Cover
4-000-522	Filter Lock
4-000-524	Filter Connector—Filter/Tubing Connector
4-000-051	Self-locking Filter, 5/pk
4-000-002	Rubber Insert, 2/pk
4-000-526	XP2 Removable Stand
4-000-036	Replacement Ni-MH Battery



### Specifications

Size	5.5"
Weight	7.8 oz.
Capacity	Accepts 1 to 100 mL pipets
Battery	Ni-MH (3.6V/700 mAh)

**Most csak 79 000 Ft+ÁFA!**



**DRUMMOND**  
SCIENTIFIC COMPANY

500 Parkway, Box 700  
Broomall, PA 19008  
Call Toll Free: 1-800/523-7480  
1-610/353-0200  
FAX: 1-610/353-6204  
Web site: [www.drummondsci.com](http://www.drummondsci.com)

6726 SZEGED TEMESVÁRI KRT 62.  
TEL: 62 / 599 750  
FAX: 62 / 599 751  
WEB: [HTTP://WWW.BIOCENTER.HU](http://www.biocenter.hu)

