

The Multi-Task Genie™ Family

THE ORIGINAL



NEW!
DIGITAL

Scientific Industries' reputation for high quality and reliability is world renowned – to date, more than 1,000,000 genuine Vortex-Genie® mixers have been used in any type of laboratory! We continue our tradition as the “Gold Standard” in vortexing and mixing with our complete Multi-Task Genie family of vortex mixers.

Vortex-Genie® 1	Vortex-Genie® 2	Vortex-Genie® 2T	Digital Vortex-Genie® 2
<ul style="list-style-type: none"> ● Violent single high-speed vortexing ● Easy-to-use Touch On operation ● Large, 6mm orbit for aggressive vortexing of difficult samples ● May be used with Pop-off Cup (shown) or optional 3-inch Platform ● Vortex-Genie quality at an even more economical price! 	<ul style="list-style-type: none"> ● OUR MOST POPULAR MODEL ● Variable speed for precise mixing from gentle to vigorous ● Broad range of attachments for most tubes, plates and other containers ● Hands-Free or Touch On control ● May be used in cold rooms and incubators 	<ul style="list-style-type: none"> ● Same great features as the Vortex-Genie 2 with an integrated timer ● Timer function: <ul style="list-style-type: none"> - Touch On: 1-60 sec. - Hands-Free: 1-60 min. or continuous ● Unit automatically stops at the end of user-programmed time ● May be used in cold rooms and incubators 	<ul style="list-style-type: none"> ● Same great features as the Vortex-Genie 2 with digital control and display of speed and time ● Digital display provides accuracy, reproducibility and repeatability ● Timer function: <ul style="list-style-type: none"> - Touch On: 1-99 sec. - Hands-Free: 1-99 min. or continuous ● May be used in cold rooms and incubators

SPECIFICATIONS

Vortex-Genie Model	Vortex-Genie 1	Vortex-Genie 2	Vortex-Genie 2T	Digital Vortex-Genie 2
Supplied	Pop-off Cup	3-inch Platform, Pop-off Cup	3-inch Platform, Pop-off Cup	3-inch Platform, Pop-off Cup
Speed (RPM)	3200 (2700 for 50Hz models)	600-3200 (2700 for 50Hz models)	600-3200 (2700 for 50Hz models)	500-3000 (2850 for 50Hz models)
Timer - Touch On: Hands-Free:	N/A	N/A	1-60 seconds 1-60 minutes or continuous	1-99 seconds 1-99 minutes or continuous
Weight	4 Kg (8.8 lb)	4 Kg (8.8 lb)	4 Kg (8.8 lb)	4 Kg (8.8 lb)
Dimensions (base) (D x W x H)	165 x 122 x 165mm (6.5 x 4.8 x 6.5in)	165 x 122 x 165mm (6.5 x 4.8 x 6.5in)	165 x 122 x 165mm (6.5 x 4.8 x 6.5in)	165 x 122 x 165mm (6.5 x 4.8 x 6.5in)
Warranty	2 years	2 years	2 years	2 years

CATALOG NUMBERS

120V, 60 Hz (North America)	SI-0136	SI-0236 (Model G560)	SI-T236	SI-A236
230V, 50 Hz - No Plug	SI-0146	SI-0246 (Model G560E)	SI-T246	SI-A246
230V, 50 Hz - European Plug	SI-0156	SI-0256 (Model G560E)	SI-T256	SI-A256
230V, 50 Hz - British Plug	SI-0166	SI-0266 (Model G560E)	SI-T266	SI-A266
230V, 50 Hz - Swiss Plug	SI-0176	SI-0276 (Model G560E)	SI-T276	SI-A276
100V, 50/60Hz	SI-0186	SI-0286	SI-T286	SI-A286

See back page for additional Vortex-Genie 2, Vortex-Genie 2T and NEW Digital Vortex-Genie 2 accessories.



ETL Listed Conforms to UL Std. 3101-1
ETL Listed Certified to CAN/CSA C22.2

U.S. Pats. 4,781,487 / 5,707,861



Scientific Industries, Inc.
Home of the Multi-Task Genie™ Family
70 Orville Drive, Bohemia, NY 11716 U.S.A.
Phone: (631) 567-4700 • Fax: (631) 567-5896
Toll Free: 888-850-6208
Customer Service: custsvc@scientificindustries.com
www.scientificindustries.com

Multi-Task Genie™ Family Accessories

A Vast Range Of Accessories For The Ultimate Multi-Task Versatility In Any Application!



Accessories for:
Vortex-Genie 2
Vortex-Genie 2T
Digital Vortex-Genie 2

Fig. 1
Model H301
Multiple Sample
Starter Set




Fig. 4
Model H302
Large Sample Set



Fig.	Catalog No.	Description	Fig.	Catalog No.	Description
1	504-0233-00	Model H301 Multiple Sample Starter Set	4	503-0282-00	Model H302 Large Sample Set
--	146-6005-00	6-inch Platform (not shown separately)	5	504-0039-00	Recessed Platform
2	504-0234-00	Microtube Foam Insert (2), 6-inch Platform Required	6	503-0278-00	9-16mm Tube Foam Inserts (2)**
3	504-0235-00	Microplate Foam Insert (2), 6-inch Platform Required	7	503-0279-00	14-29mm Tube Foam Inserts (2)**
11	146-3011-00	Pop-off Cup	8	503-0280-00	29-37mm Tube Foam Inserts (2)**
12	0K-0500-902	3-inch Platform with Rubber Cover	9	503-0277-00	Blank Foam Inserts (2)**
13	SI-0563	TurboMix™ Attachment, 1.5ml Tubes	--	503-0283-00	Assorted Foam Inserts (4)**
--	SI-0562	TurboMix™ Attachment, 2.0ml Tubes (not shown)	10	568-0001-00	Elastic Bands (2)**
--	0A-0563-010	1.5ml Microtube Holder* (not shown separately)	--	SI-0510	Two-Tier Microplate Foam Insert (not shown)
14	0A-0563-011	2.0ml Screw-cap Microtube Holder*	15	SI-0511	Large Ampule/Tube Attachment
			16	SI-0570	Small Ampule/Tube Attachment

* May be used with TurboMix attachment or Disruptor Genie™

** Recessed Platform Required

NOTE: Vortex-Genie 1 may only be used with Pop-off Cup and 3-inch Platform

U.S. Pat. 5,707,861



Scientific Industries, Inc.
Home of the Multi-Task Genie™ Family
70 Orville Drive, Bohemia, NY 11716 U.S.A.
Phone: (631) 567-4700 • Fax: (631) 567-5896
Toll Free: 888-850-6208
Customer Service: custsvc@scientificindustries.com
www.scientificindustries.com