

JB Nova, General purpose unstirred water baths

DESCRIPTION

The reliability, quality and consistent performance of Grant products have made Grant a leading manufacturer of water baths for decades.

- Set and Forget™ technology - fast heat-up, reliable temperature control
- Simple to use controls, clear bright display
- 3 year warranty. UK design and manufacture
- Display lock - disables front panel controls
- Drain tap (not on 5L version), tray and lid included
- Advanced dry start and run dry protection - prevents costly service repairs



SPECIFICATIONS

Tank capacity	5 / 12 / 18 / 26 L
Temperature range	ambient +5 to 95°C
Temp. display and setting resolution	0.5°C
Temp stability (DIN 12876) @ 70°C	±0.5°C
User defined calibration points	1
Nominal operating voltage	120 or 230 V

CAT. NUMBER

JBN5	230VAC 50/60Hz Euro plug
JBN5 US	120VAC 50/60Hz US plug
JBN12	230VAC 50/60Hz Euro plug
JBN12 US	120VAC 50/60Hz US plug
JBN18	230VAC 50/60Hz Euro plug
JBN18 US	120VAC 50/60Hz US plug
JBN26	230VAC 50/60Hz Euro plug
JBN26 US	120VAC 50/60Hz US plug